

# Tecella Pico

## USB-powered Miniature Patch Clamp Amplifier

Pico is a feature-rich, low-noise patch clamp amplifier with an integrated digitizer and headstage. Pico is ideally suited for whole cell (Vclamp and Iclamp) and single-channel (patch, planar lipid bilayer and synthetic nanopore) recording, as well as cellular electrochemistry.



Clean Head Switching™ technology allows software controlled switching between voltage clamp and current clamp, or between multiple feedback resistors, without introducing any artifacts. This enables one Pico to support a wide range of applications including Whole Cell, Single Channel, Multi-Cell, Bilayers, Electrochemistry, and Current Clamp.

In voltage clamp mode, the Pico provides 5 feedback gain resistors ranging from 10 M $\Omega$  to 10 G $\Omega$ . In current clamp mode, the Pico provides 3 range settings from  $\pm 2$  nA to  $\pm 200$  nA.

Internal Model Cell allows for self calibration of the experiment setup, as well as, rapid post-experiment amplifier validation.




Software support for the Pico includes TecellaLab acquisition software and WinWCP software. Additional 3<sup>rd</sup> party software support is planned in the near future. SDK / API is provided for custom development and system integration.

RMS noise is 1 pA at 1 G $\Omega$  gain setting.

- Fully featured 1-channel amplifier
- USB powered
- Voltage Clamp and Current Clamp
- Series resistance compensation
- Multiple capacitance compensations
- Very small size
- Integrated headstage and digitizer
- Under 1 Watt power consumption
- Low Noise

The Patented Spread Frequency Compensation can automatically compensate any arbitrary capacitance profile in approximately 3 seconds.



USB cable is all that is needed to power-up, control, and record from the Pico.

One BNC-type electrode holder is included with the Pico. Any BNC-type electrode holder can be used with the Pico.



Sutter Instruments, Siskiyou, and NeoBiosystems provide hardware to mount the Pico directly onto their micro manipulators.

## Accessories



Included

USB Cable



Electrode Holder



Break-Out Box  
for Triggers and  
External DAQ Interface

## Specifications

Integrated Digitizer	20 kHz sampling rate 16-bit A/D (18-bit internal resolution) Stimulus voltage ranges: ±250 mV ±2000 mV Zap voltage range of ±1000 mV
Feedback Gain Settings	10 MΩ, 100 MΩ, 1 GΩ, 3.3 GΩ, 10 GΩ
Low RMS Noise (DC to 3kHz)	0.3 pA @ 10 GΩ, 1.0 pA @ 1 GΩ, 7 pA @ 100 MΩ
Filters	Programmable 2-pole Low-Pass Filter (analog hardware circuit) Digital Filter available in TecellaLab software
Compensations	Up to 4 Capacitance Compensations Cfast x 1, Cslow x 3 0-100 pF per compensation Series Resistance Compensation Offset Compensation (±250 mV) Optional Active Leak Compensation
Current Clamp	±2 nA range with 1.25 pA resolution ±20 nA range with 12.5 pA resolution ±200 nA range with 125 pA resolution
Break-Out Box (optional)	For controlling external devices via the Pico's integrated digitizer: Digital Out x 4 (For Triggers)  For interfacing Pico to 3 <sup>rd</sup> Party Digitizers (DAQs): Analog In Analog Out Gain Telegraph 2 Current Clamp Control signals - "Imode" puts Pico in I=0/Iclamp mode - "Iclamp" engages Iclamp mode Switch to enable 3 <sup>rd</sup> Party Digitizer control
Computer Interface	USB
Software	WinWCP by the University of Strathclyde jClamp by SciSoft Company TecellaLab software, with Data Export to ATF, tab-formats SDK/API available Other 3 <sup>rd</sup> Party Software support planned
Mechanical & Power	USB powered. Does not require a separate power supply. 5.8 in x 1.8 in x 0.7 in (14.5 cm x 4.6 cm x 1.7 cm) 3 ounces (85 grams)
Power Consumption	0.8 Watt

U.S. Patents 7,741,829  
8,163,147  
8,163,527

97 Chaumont Cir  
Foothill Ranch, CA 92610  
USA  
Phone: +1-714-641-1709  
Fax: +1-714-641-1569  
E-Mail: info@tecella.com